

Original Research Article

Forced Degradation Behaviour with a Developed and Validated RP-UPLC Method for Estimation Lamivudine and Dolutegravir in Combined Pharmaceutical Dosage Form.

ABSTRACT:

A simple, fast, precise, specific, and accurate reversed phase Ultra Performance liquid chromatographic (UPLC) method was developed and validated for the forced degradation studies of the Lamivudine and Dolutegravir in tablet dosage form. Chromatogram was run HSS C18 (2.6 x 50 mm, 1.6 µm). Mobile phases containing 70% 0.01N disodium hydrogen phosphate: 30% Methanol). The elution of analytes was achieved with a flow rate at 0.3 mL/min. Potassium dihydrogen phosphate was used as a buffer in this experiment. Temperature was maintained at 30°C. Optimized wavelength selected at 260 nm. The detector response was linear in the concentration range of 25- 150 µg/mL respectively. Retention time of Lamivudine and Dolutegravir were found to be 1.408 min and 1.739 min. %RSD of the Lamivudine and Dolutegravir were and found to be 0.8 and 0.8 respectively. %Recovery was obtained as 100.39% and 100.37% for Lamivudine and Dolutegravir respectively. LOD, LOQ values obtained from regression equations of Lamivudine and Dolutegravir were 0.41, 1.25 and 0.09, 0.26 respectively. Regression equation of Lamivudine is $y = 24270x + 12218$ and $y = 34783x + 1060$ Dolutegravir. Since retention times and run times were reduced, the method established was easy and premium and it could be used in periodic quality control tests in industries.

KEY WORDS: Lamivudine, Dolutegravir, RP-UPLC, Forced degradation.

INTRODUCTION:

Dolutegravir is Chemically (4R,12aS)-N-[(2,4-Difluorophenyl)methyl]-3,4,6,8,12,12a-hexahydro-7-hydroxy-4-methyl-6,8-dioxo-2H-pyrido[1',2':4,5]pyrazino[2,1-b][1,3]oxazine-9-carboxamide (Figure No:1). Its molecular Formula $C_{20}H_{19}F_2N_3O_5$ and Molecular weight 419.38 g/mol¹. Dolutegravir is an HIV-1 integrase strand transfer inhibitor of the second generation. For the treatment of HIV infection, dolutegravir is now in Phase III clinical studies. Dolutegravir has been proven to inhibit HIV replication in cells infected with a self-inactivating PHIV lentiviral vector, such as peripheral blood mononuclear cells (PBMCs), MT-4 cells, and CIP4 cells. It is slightly soluble in water and Acetonitrile. Dolutegravir is an HIV-1 antiviral agent. It binds to the active site of HIV integrase and prevents the strand transfer stage of retroviral DNA integration. This is a critical phase in the HIV replication cycle that will result in viral activity being inhibited. Dolutegravir has a mean EC50 value of 0.5 nm (0.21 ng/mL) to 2.1 nm (0.85 ng/mL) in peripheral blood mononuclear cells (PBMCs) and MT-4 cells. Dolutegravir is authorised for the treatment of human immunodeficiency virus type 1 (HIV-1) infection in adults and children aged 12 years and older in conjunction with other antiretroviral medicines and weighing at least 40 kg².

Lamivudine chemically known as 4-amino-1-[(2R, 5S)-2-(hydroxymethyl)-1,3-oxathiolan-5-yl]-1,2-dihydropyrimidin-2-one. Its molecular formula $C_8H_{11}N_3O_3$ (Figure No:2). Lamivudine is a synthetic nucleoside analogue that is phosphorylated intracellularly to produce lamivudine triphosphate, an active 5'-triphosphate metabolite (L-TP). HIV reverse transcriptase and HBV polymerase integrate this nucleoside analogue into viral DNA, resulting in DNA chain termination. Lamivudine is a nucleoside reverse transcriptase inhibitor (NRTI) that inhibits viral DNA synthesis in the Human Immunodeficiency Virus Type 1 (HIV-1) and Hepatitis B (HBV). Lamivudine can produce active metabolites when phosphorylated, which compete for inclusion into viral DNA. Lamivudine metabolites function as a chain terminator of DNA synthesis by competitively inhibiting the activity of the HIV reverse transcriptase enzyme via DNA incorporation. Incorporated nucleoside analogues prohibit the development of a 5'-OH group because they lack a 3'-OH group. Phosphodiester linkage that is essential for DNA chain elongation²⁻⁴.

Extensive literature survey was carried out which revealed that there is no work carried out especially on Forced degradation method of the Lamivudine and Dolutegravir in tablet dosage form using UPLC method. The specific aim of the research was to develop a UPLC

method and performing forced degradation studies of Lamivudine and Dolutegravir in bulk and formulated dosage form and to validate the proposed methods in accordance with ICH guidelines for the intended analytical application⁵⁻¹²

MATERIALS AND METHODS:

Instruments and Apparatus:

Standard Lamivudine and Dolutegravir pure drugs (API), Combination Lamivudine and Dolutegravir tablet (Dovato), Distilled water, Acetonitrile, Phosphate buffer, Methanol, Potassium dehydrogenate ortho phosphate buffer, Ortho-phosphoric acid. All the above chemicals and solvents are from Rankem Pvt. Ltd, India.

A Ultra-High Performance liquid chromatography (WATERS UPLC Auto Sampler TUV Detector with Empower 2 Software). UV-VIS spectrophotometer PG Instruments T60 with special bandwidth of 2 mm and 10 mm and matched quartz cells integrated with UV win 6 Software was used for measuring absorbances of Lamivudine and Dolutegravir solutions.

Preparation of Standard stock solutions:

Accurately weighed 75 mg of Lamivudine, 12.5 mg of Dolutegravir and transferred to individual 50 mL volumetric flasks separately. Both flasks were sonicated for 10 minutes after 3/4 of the diluents were introduced. Standard stock solution 1 and 2 were prepared in flasks with diluents and labelled. (1500 µg/mL of Lamivudine and 250 µg/mL of Dolutegravir).

Comment [A1]: Maintain accurate space in between two specific words and digits.

Correct this in all remain article.

Preparation of Standard working solutions (100% solution):

1mL from each stock solution was Pipetted out and taken into a 10 mL volumetric flask and made up with diluent. (150 µg/mL Lamivudine of and 25 µg/mL of Dolutegravir).

Preparation of Sample stock solutions:

10 tablets were weighed and was transferred into a 100 mL volumetric flask, 50 mL of diluents was added and sonicated for 25 min, further the volume was made up with diluent and filtered by HPLC filters (3000 µg/mL of Lamivudine and 500 µg/mL of Dolutegravir).

Preparation of Sample working solutions (100% solution):

0.5 mL of filtered sample stock solution was transferred to 10mL volumetric flask and made up with diluent (150 µg/mL of Lamivudine and 25 µg/mL of Dolutegravir).

Buffer preparation:

Accurately weighed 1.41gm of sodium dihydrogen Ortho phosphate in a 1000 mL of Volumetric flask add about 900 mL of milli-Q water added and degas to sonicate and finally make up the volume with water.

Method Development:

Chromatographic analysis was performed on Acquity UPLC HSS C18 (2.6 x 50 mm, 1.6 µm). The mobile phase consists of 70% 0.01N Na₂HPO₄: 30% Methanol was used through the analysis. The flow rate was 0.3 mL/min, the injection volume was 1.0mL, column temperature was 10°C and detection was performed at 260 nm using a UV detector.

METHOD VALIDATION:

The analytical method was verified for linearity, accuracy, precision, specificity, robustness, and ruggedness in accordance with ICH recommendations.

Linearity:

Linearity was demonstrated from 1500 µg/mL of Lamivudine and 250 µg/mL of Dolutegravir of standard concentration using minimum six calibration levels (25%, 50%, 75%, 100%, 125% and 150%) for both the drugs. The data was analysed using the linear regression approach.

Accuracy:

The proximity between the reference value and the obtained value is expressed by the accuracy of an analytical method. The accuracy of the method was evaluated in triplicate at three concentration levels, 50%, 100%, and 150% of standard solutions of Lamivudine and Dolutegravir.

Robustness:

Robustness conditions like Flow minus (0.2 mL/min), Flow plus (0.4 mL/min), mobile phase minus, mobile phase plus, temperature minus (25°C) and temperature Plus (35°C) was

Comment [A2]: Use only one fullstop after completeing sentence.

maintained and samples were injected in duplicate manner. The system suitability parameters were relatively unaffected, and all of them were passed. %RSD was within the limit.

Specificity:

Checking of the interference in the optimized method. We should not find interfering peaks in blank and placebo at retention times of these drugs in this method. So this method was said to be specific and excellent.

Precision:

Preparation of Sample stock solutions:

10 tablets were weighed and was transferred into a 100 mL volumetric flask, 50 mL of diluents was added and sonicated for 25 min, further the volume was made up with diluent and filtered by HPLC filters (3000 µg/mL of Lamivudine and 500 µg/mL of Dolutegravir).

Preparation of Sample working solutions (100% solution):

In a 10 mL volumetric flask, 0.5 mL of filtered sample stock solution was transferred and made up with diluent and made up with diluent. (150 µg/mL of Lamivudine and 25 µg/mL of Dolutegravir).

Forced degradation studies:

Oxidation:

To 1 mL of stock solution of Lamivudine and Dolutegravir, 1 mL of 20% hydrogen peroxide (H₂O₂) was added separately. The solutions were maintained at 60°C for 30 minutes. For UPLC study, the resultant solution was diluted to obtain 150 µg/mL & 25 µg/mL solution and 1.0 µL were injected into the system.

Acid Degradation Studies:

To 1 mL of stocks solution Lamivudine and Dolutegravir, 1 mL of 2N Hydrochloric acid was added and refluxed for 30 mins at 60⁰c. The resultant solution was diluted to obtain 150 µg/mL & 25 µg/mL solution and 1.0 µl solutions were injected into the system.

Alkali Degradation Studies:

To 1 mL of stock solution Lamivudine and Dolutegravir, 1 mL of 2N sodium hydroxide was added and refluxed for 30mins at 60⁰c. The resultant solution was diluted to obtain 150 µg/mL & 25 µg/mL solution and 1.0 µL were injected into the system.

Dry Heat Degradation Studies:

To investigate dry heat degradation, the standard medication solution was baked for 6 hours at 105°C. The resulting solution was diluted to 150 µg/mL & 25 µg/mL for UPLC testing. To and 1.0 µL were injected into the system.

Photo Stability studies:

The photochemical stability of the drug was also studied by exposing the 1500 µg/mL & 250 µg/mL solution to UV Light by keeping the beaker in UV Chamber for 7 days or 200-Watt hours/m² in photo stability chamber. For UPLC study, the resultant solution was diluted to obtain 150 µg/mL & 25 µg/mL solutions and 1.0 µL were injected into the system.

Neutral Degradation Studies:

Stress testing under neutral conditions was studied by refluxing the drug in water for 6 hrs at temperature of 60°. For UPLC study, the resultant solution was diluted to 150 µg/mL & 25 µg/mL solution and 1.0 µL were injected into the system. (Figure 6 to 9), (Table 6)

RESULTS AND DISCUSSION:**Optimized method:**

Lamivudine and Dolutegravir were eluted at 1.408 min and 1.739 min respectively with good resolution. Because the plate count and tailing factor were both good, this method was refined and validated. (Figure No: 3)

Linearity:

Six linear concentrations of Lamivudine (37.5-225 µg/mL) and Dolutegravir (6.25-37.5 µg/mL) were injected in a duplicate manner. Average areas were mentioned above and linearity equations obtained for Lamivudine was $y = 24270x + 12218$ and of Dolutegravir was $y = 34783x + 1060$. The two drugs have a correlation value of 0.999. (Figure No: 4 & 5), (Table 1).

Precision:

From a single volumetric flask of working standard solution six injections were given and the obtained areas were mentioned above. Average area, standard deviation and % RSD were calculated for two drugs. % RSD obtained as 0.8% and 0.8% respectively for Lamivudine and Dolutegravir. (Table 5).

Accuracy:

Three levels of Accuracy samples were prepared by standard addition method. Triplicate injections were given for each level of accuracy and mean %Recovery was obtained as 100.39% and 100.37% for Lamivudine and Dolutegravir respectively.(Table 2 and 3)

Robustness:

Robustness conditions like Flow minus (0.2 mL/min), Flow plus (0.4 mL/min), mobile phase minus (75B:25M), mobile phase plus (65B:35M), temperature minus (25°C) and temperature plus (35°C) was maintained and samples were injected in duplicate manner. The system suitability parameters were largely unaffected, and all of them were passed. %RSD was within the limit. (Table 4).

CONCLUSION

Comment [A3]: Write conclusion specific and detail.

A simple, Accurate, precise method was developed for the simultaneous estimation of the Lamivudine and Dolutegravir in Tablet dosage form. Retention time of Lamivudine and Dolutegravir were found to be 1.408 min and 1.739 min. %RSD of the Lamivudine and Dolutegravir were and found to be 0.8 and 0.8 respectively. %Recovery was obtained as 100.39% and 100.37% for Lamivudine and Dolutegravir respectively. LOD, LOQ values obtained from regression equations of Lamivudine and Dolutegravir were 0.41, 1.25 and 0.09, 0.26 respectively. Regression equation of Lamivudine is $y = 24270x + 12218$ and $y = 34783x + 1060$ of Dolutegravir. The technique adopted was convenient and cost-effective, and it could be used in routine quality control tests in industries.

COMPETING INTERESTS DISCLAIMER:

Authors have declared that no competing interests exist. The products used for this research are commonly and predominantly use products in our area of research and country. There is absolutely no conflict of interest between the authors and producers of the products because we do not intend to use these products as an avenue for any litigation but for the advancement of knowledge. Also, the research was not funded by the producing company rather it was funded by personal efforts of the authors.

REFERENCE:

1. Indian Pharmacopoeia, Indian Pharmacopoeial Commission (IPC), Ministry of Health and Family Welfare: Ghaziabad, Delhi.2014; Vol. 1: 1900-1902.
2. A to Z drug facts. Ed. 7; Wolters Kluwer Health Inc. 2006; 7(2):1313-1387.

3. Tripathi K. D, Essentials of Medical Pharmacology. JP Medical Ltd: 2006; 1(6): 565-568.
4. ICH Guidelines: Text on validation of analytical procedures. In proceedings of International Conference on Harmonization.1994; 1-5.
5. Somshankar D, Mahesh D. Simultaneous estimation of Lamivudine, Abacavir and Dolutegravir by UPLC method. International journal of applied pharmaceutics. 2018; 1(10): 44-49. doi.org/10.22159/ijap.2018v10i1.21156
6. Devanna N, Benzil D, Ramachandraiah C. Development and validation for the simultaneous estimation of Dolutegravir and Lamivudine in drug product by RP-HPLC. International journal of research in pharmaceutical and nano sciences. 2017; 6(4):173 - 180.
7. Michael ES, UPLC: An Introduction and Review Waters Corporation, Milford, Massachusetts, USA. Journal of Liquid Chromatography & Related Technologies. 2005; 28(3): 1253–1263.
8. Lars Y, Honore HS. On-line turbulent-flow chromatography high-performance liquid chromatography mass spectrometry for fast sample preparation and quantitation. Journal of Chromatography. 2003; 1(10):59–67.
9. Goodwin L, White SA, Spooner N. Evaluation of ultra-performance liquid chromatography in the bioanalysis of small molecule drug candidates in plasma. Journal of Chromatography Science. 2007; 45(6): 298–304.
10. Sowjanya B, Development and validation for the simultaneous estimation of Lamivudine and Dolutegravir in drug product by RP-HPLC, indo American journal of pharmaceutical research, 2018; 12(3):225-229.
11. Dhanve A, Rao JR, Dhale C. Development And Validation Of Stability Indicating HPTLC Method For Simultaneous Estimation Of Lamivudine And Dolutegravir Sodium In Bulk And Pharmaceutical Dosage Formulation. International Journal Of Pharmaceutical Sciences And Research. 2018; 9(11):4701-8.
12. Kalpana N, Shanmukha KV, Ramachandran D. Development and Validation for the Simultaneous Estimation of Lamivudine, Tenofovir Disoproxil and Dolutegravir in Drug Product by RP-HPLC. Journal of Pharmaceutical Science and Research. 2017;

Comment [A4]: Don't use extra space for reference. From ref. no. 11 to 12 its having huge distance. Kindly maintain proper space.

FIGURES:

Comment [A5]: It's better to attach tables and figures are as per specific points. It will easy to understand.




Figure 1.Structure of Dolutegravir




Figure 2.Structure of Lamivudine




Figure 3. Optimized Chromatogram




Figure 4. Calibration curve of Lamivudine




Figure 5. Calibration curve of Dolutegravir




Figure 6. Acid degradation chromatogram




Figure 7. Base degradation chromatogram




Figure 8. Peroxide degradation chromatogram




Figure 9. Thermal degradation chromatogram

TABLES:

Table 1. Linearity table for Lamivudine and Dolutegravir

Lamivudine		Dolutegravir	
Conc (µg/mL)	Peak area	Conc (µg/mL)	Peak area
0	0	0	0
25	922806	25	210556
50	1878662	50	445258
75	2712498	75	655630

100	3638097	100	872691
125	4569129	125	1084994
150	5477284	150	1303577

Table 2. Accuracy of Lamivudine

% Level	Amount Spiked (µg/mL)	Amount recovered (µg/mL)	% Recovery	Mean %Recovery
50%	75	74.63	99.50	100.39%
	75	74.93	99.91	
	75	75.55	100.74	
100%	150	147.86	98.58	
	150	151.98	101.32	
	150	150.99	100.66	
150%	225	228.58	101.59	
	225	226.97	100.88	
	225	225.78	100.34	

Table 3. Accuracy of Dolutegravir

% Level	Amount Spiked (µg/mL)	Amount recovered (µg/mL)	% Recovery	Mean %Recovery
50%	12.5	12.49	99.94	100.37%
	12.5	12.48	99.83	
	12.5	12.54	100.34	
100%	25	25.34	101.36	
	25	25.28	101.12	
	25	25.25	101.00	
150%	37.5	37.32	99.51	
	37.5	37.69	100.50	
	37.5	37.38	99.68	

Table 4. Robustness data for Lamivudine and Dolutegravir

S.No	Condition	%RSD of Lamivudine	%RSD of Dolutegravir
1	Flow rate (-) 0.2mL/min	0.7	1.0
2	Flow rate (+) 0.4mL/min	0.3	0.5
3	Mobile phase (-)	0.2	0.4

	75B:25M		
4	Mobile phase (+) 65B:35M	0.7	1.1
5	Temperature (-) 25°C	0.1	0.1
6	Temperature (+) 35°C	0.2	0.7

Table 5. System precision table of Lamivudine and Dolutegravir

S. No	Area of Lamivudine	Area of Dolutegravir
1.	3640080	882110
2.	3662158	882532
3.	3623764	870302
4.	3646966	872703
5.	3650353	878651
6.	3710141	889046
Mean	3655577	879224
S.D	29577.7	6902.1
%RSD	0.8	0.8

Table 6. Degradation data for Lamivudine and Dolutegravir

Type of degradation	Lamivudine			Dolutegravir		
	Area	Percentage Recovered	Percentage Degraded	Area	Percentage Recovered	Percentage Degraded
Acid	3408992	93.16	6.84	832130	94.55	5.45
Base	3496767	95.56	4.44	832826	94.63	5.37
Peroxide	3520360	96.20	3.80	814105	92.50	7.50
Thermal	3623820	99.03	0.97	846632	96.20	3.80
UV	3617941	98.87	1.13	849805	96.56	3.44

Water	3651379	99.79	0.21	875607	99.49	0.51
-------	---------	-------	------	--------	-------	------

UNDER PEER REVIEW